

White-Light Scanning Digital Microscopy

J. Paul Robinson

Purdue University, West Lafayette, Indiana

Ben Gravely

Cosmic Technologies Corporation, Raleigh, North Carolina

INTRODUCTION

Digital microscopy has for many years been a technology difficult to implement. Finally, however, it is clear that new technologies are beginning to reshape the field of microscopy. The dramatic increase in the performance and cost effectiveness of computers, digital storage systems, and telecommunications is fueling rapid growth in the creation and storage of digital microscope images. There are several important issues that must be addressed in dealing with digital images: (1) the image collection system, which is usually a charge-coupled device (CCD) camera; (2) the image capture computer, including necessary boards and software; and (3) image manipulation (image processing) and data transfer from one site to another.

The development of digital microscopy is being driven by the need to make critical decisions based upon the rapid collection and analysis of images. For fifty years a desire to use imaging to define cancerous cells from normal has been a motivating force (Mellors et al., 1952). It has been apparent for many years that these decisions should be made with the aid of imaging technologies in an automated fashion (Bartels et al., 1984). Early attempts at using rapid imaging technologies

(video) proved very difficult and were unable to produce high-resolution images (Bibbo et al., 1983; Ingram and Preston, 1964; Prewitt and Mendelsohn, 1966) or simply reached an impasse because of lack of rapid imaging technology (Ingram and Preston, 1970). With the development of more powerful computers and automated microscopes, complete systems for diagnostic use began to take shape (Ploem et al., 1986, 1988). The trend today in medicine, for example, is the move toward electronic data integration of all patient records, including diagnostic images, and this promises to revolutionize not only diagnosis but also teaching with archival databases of images, automated analysis, publication, and local and remote consultation.

Digital microscopy is a key element in a combination of technologies that has made possible the development of telepathology, which involves the electronic transfer of pathological images between geographically distant locations for medical consultation. The remote control of a distant microscope and the transfer of high-quality images over normal telecommunications systems coupled with text and voice will have a significant impact on clinical pathology, resulting in more efficient diagnosis and a greater availability of quality medical services to remote areas.

Telepathology eliminates the need to physically transport critical samples through delivery services, whose availability and reliability may be nonexistent in many areas of the world. Telepathology can be divided into two forms. Still telepathology involves sending "snapshots" of areas of interest from the client to the consultant and works best in a consultation environment between peers. The client should be qualified to select the fields of interest. This technology frequently utilizes high-quality image capture cameras, and images are downloaded in near real time (sometimes several minutes) for remote evaluation, such as that recently reported by Singson et al. (1999). Other less advanced systems have been proposed using regular telephone lines and modems; however, these suffer from loss of high-resolution images and significant and unacceptable time delays (Vazir et al., 1998). Dynamic telepathology involves the robotic operation of a remote microscope by the consulting pathologist, who selects the fields and focuses just as he or she would with a local microscope. A number of publications have appeared recently demonstrating remote microscope control using the internet as a somewhat universal platform (Nagata and Mizushima, 1998; Szymas and Wolf, 1998). The remote microscope responds to the commands and relays images back to the pathologist. Only a qualified technician is required at the remote microscope to prepare slides and place them in the system. Dynamic telepathology is the most demanding of all digital pathology tasks because the high image quality and speed of operation required are very difficult to achieve. Before dynamic telepathology can become common, many complex issues must be solved in image capture and transmission systems.

This chapter discusses one potential solution to the technology problems—a white-light flying-spot scanning microscope system. This microscope is an integrated, fully digital color microscope that utilizes a flying-spot light source, advanced image-processing techniques, and simple straightforward operation. Although the fundamental principles upon which this technology is based were developed nearly a

half century ago, only now have technological advancements allowed it to become an economic reality.

SCANNING SPOT

The color scanning microscope (COSMIC) is based on spot-scanning techniques developed in the 1950s for various applications, including film-to-video conversion for television and even black-and-white video microscopes. Figure 13.1 shows the optical diagram of the instrument. A specially designed cathode ray tube (CRT) has a




Fig. 13.1. Schematic layout of the COSMIC system showing the light source, optical path, and signal collection components. The spot scans across the CRT create a point of light that illuminates only the Airy disk area on the sample. This reduces bleaching and scattered light on the sample. Since the illumination point is the smallest that can be created by the objective, the system resolution is limited only by the optical resolution of the objective, not the detector components.

white raster pattern on its screen. This pattern is imaged through a microscope relay lens and objective lens onto the sample placed on the microscope stage.

The size of the raster image on the sample is determined by the lens system magnification. For example, an ordinary microscope may have a 40-power objective lens and a 10-power eyepiece lens, giving a combined magnification of 400 from the sample to the eye. Operated in reverse, an object placed outside the eyepiece lens will be demagnified 400 times down to the sample plane.

This one-to-one mapping of the scanning spot from the CRT screen onto the sample plane creates a tiny white probe in a dark field that scans over the sample. The spot characteristics on the sample are completely determined by the microscope objective. As a result, the probing spot has the smallest diffraction-limited diameter possible for each objective lens.

The light transmitted through the sample is divided by dichroic mirrors into red (580–650 nm), green (520–580 nm), and blue (450–520 nm) spectral components and passed to three detectors. A video display reconstructs the image of the sample using the x - y position of the scanning spot and the three color intensities. In the simplest case, the scan frequency of the CRT raster matches the scan frequency of the display, giving a direct correspondence between the spot position on the CRT and the spot position on the display.

DETECTION SYSTEM

The detectors in this digital microscope are not the CCD imaging devices used in most “digital” cameras today. Instead, the detectors are three nonimaging photomultiplier tubes (PMTs) that collect the intensity signals simultaneously in the three spectral bands. Color registration is inherently perfect since there is only one spot scanning over the sample and the PMT detectors read intensity only in the color bands. The PMTs used are ordinary types available from several manufacturers, capable of achieving full color saturation with signals of less than 100 nW of power. More expensive or cooled PMT detectors could be used for special applications.

SIGNAL PROCESSING

Performance can be improved by reducing the scan frequency of the CRT well below the scan frequency of the display. The data are read in at one rate, then presented to the display at a faster rate. The choice of scan rate at the sample is determined by the number of pixels in the image and the desired optical contrast function. From a purely technical standpoint, the number of pixels in the digital image could be chosen as any value, but practical considerations of the data file size, availability of displays, processing and transmission times, and correlation with direct viewing set reasonable boundaries to the image file size.




Fig. 13.2. Signal processing path showing how the transmitted light is handled for each spectral channel.

These considerations led to the choice of image size for COSMIC as 1280×1024 pixels with 24 bits/pixel color depth. Displays are readily available for this resolution and the data set is manageable for computation, storage, and transmission. The image size is a reasonable match to the number of optical resolution elements in various microscope images as viewed directly by the eye. The selected frame rate of 13.4 Hz is fast enough to observe live cells. The resulting data rate is 420 Mbits/s, which is at the high end of current imaging technologies.


Figure 13.2 illustrates the electronic signal path. The three color signals from the detectors are amplified and sent to a digital imaging board that performs multiple functions. First, the three signals are digitized into parallel 1280×1024 math buffers. The data in the three math buffers are transferred to three display buffers, converted into analog signals, and sent to the display at a 60-Hz frame rate. Math

functions for averaging up to 256 images and summing (integrating) up to 9999 images are executed in the math buffers by programming built into the imaging board. These functions operate at the full frame rate.


ZOOM MAGNIFICATION

The optical magnification can be changed instantaneously up to 300% ($3\times$) by changing the scanned area on the sample. The principle is illustrated in Figure 13.3. The scan raster overlaying the sample is shown on the left and the resulting image on the display is shown on the right. The image is always digitized into 1280×1024 pixels, so the data set remains constant for any magnification, and the pixel size never changes within the field of view.

Operating an objective lens at three times its rated magnification will exceed its classical resolution limit, resulting in so-called empty magnification, which has tra-



(a) Unzoomed Scan



(b) Zoomed Scan

Fig. 13.3. As the scan area decreases, the specimen fills more of the image window. The display area is constant in size, thus the specimen is enlarged. The image is always digitized at 1280×1024 data points, creating a constant file size at all zoom values.



Fig. 13.4. Example of a zoomed image. Here, a mosquito head is imaged using the microscope's 10x objective with no zoom (A), zoom at 1.5x (B), and zoom at 3x (C). *See color plates.*

ditionally been considered a useless exercise. However, what really happens is that each optical resolution element may be digitized into as many as six data points, resulting in a process called digital oversampling. Oversampling tends to push resolution limits beyond standard conventions and results in "superresolution." Oversampling is discussed in more detail below.

The zoom feature also reduces the need for larger digital image files having more pixels since one can instantaneously zoom in to display maximum resolution at reduced field size. COSMIC's constant file size means that transmission times are invariant for all images. Figure 13.4 shows an example of a zoomed image. Here, a mosquito head is imaged using COSMIC's 10x objective with no zoom (A), zoom at 1.5x (B), and zoom at 3x (C).

BRIGHTNESS AND FOCUS CHANGE

Some unexpected characteristics occur with the zoom capability. In an ordinary optical zoom lens two things happen when the field size is changed. The brightness changes because the numerical aperture (NA) has changed (the change in focal length changes the NA), and the focus changes because it is virtually impossible to make a variable-focus lens with moving parts that can maintain constant focus throughout its range. In a point scanning microscope, on the other hand, the zoom feature does not change the focus or the brightness, which are not simply constants but universal constants. The focus does not change because mechanical distances are fixed. There are no alterations in the distance between the lenses and image planes at any time. The changes are related only to the area scanned on the sample. The size of the scanned area has been decreased, but all the mechanical distances remain the same. Another important characteristic of the zoom feature is that the brightness does not change. The illumination energy density on the sample varies with the size of the scanned area. For a large area, the energy per unit area is low. This is because the illumination emitted by the CRT is constant. Only the scanned area changes and therefore the energy density changes linearly with area.

The change in brightness between two different illumination areas can be shown mathematically to be inversely proportional to the ratio of the areas. Since the illumi-

nation on the sample changes proportionally to the area ratio, the brightness in the final image does not change with zoom.

CONTRAST ENHANCEMENT AND SUPPRESSION

In other spot-scanning instruments, such as laser scanning microscopes and electron beam microscopes, the intensity of the spot is constant as it scans over the sample. With COSMIC, however, it is possible to modulate the brightness of the spot for every pixel in the sample plane. If the spot is modulated using data from the sample, then several interesting effects are possible. Figure 13.5A illustrates the usual condition of a constant-intensity scan overlaid on the sample and the resultant image on the display. In Figure 13.5B, the CRT raster has a positive image of the sample superimposed on it. When transferred to the sample plane, the spot goes bright when the sample is bright




Fig. 13.5. Concepts of spot modulation demonstrated. The brightness of the scanning spot is varied point by point as it traces over the sample in order to increase or decrease the light/dark signal swing. Three effects are shown: standard (A), contrast enhancement (B), and contrast suppression (C).



Fig. 13.6. Image enhancement features of COSMIC. The normal-mode image of Figure 4 (A), the image in contrast suppression mode (B), and the image in color-inverted mode (C). This mode is particularly useful to identify thin structures in cells where the color inversion highlights the objects. *See color plates.*

and dark when the sample is dark, resulting in increased contrast in the image on the display. If the CRT has a negative image superimposed on it (Figure 13.5C), then the spot is dark when the sample is bright and bright when the sample is dark, suppressing contrast in the image. The contrast enhancement feature is very useful for examining unstained samples, and contrast suppression is useful for overstained or dark samples. Many of these samples cannot be successfully imaged with conventional microscopes. Figure 13.6 (A) is the same normal-mode sample as Figure 4. Figure 13.6 (B) demonstrates the contrast enhancement due to spot modulation, and Figure 13.6 (C) shows digital color inversion.

RESPONSE FUNCTION

Most television imaging and display tubes have a nonlinear response represented by an exponential equation using the variable “gamma.” Gamma correction circuits are used routinely in video amplifiers to linearize the response.

Although the response function of the digital microscope is inherently linear, a programmable response function has been designed into the imaging board. In terms of gray-level histograms, spot modulation expands or compresses the whole histogram, while the gamma can adjust the midlevel values and not the end points. The gamma response function and spot modulation functions can be used together to fine tune the desired image enhancement effects.

In addition to the modes discussed, digital brightness, contrast, and background correction can be applied. The presence of these functions recognizes the need for simple, real-time image processing, which reduces the need for postprocessing and ensures that the image saved is of the highest quality.

AVERAGING AND SUMMING

In order to boost weak and smooth noisy signals, two digital processing modes have been implemented in hardware using the math buffer. The averaging mode performs a running average of up to 256 images in steps of 2^n , where $n = 1, \dots, 8$. It is also possible to sum (or integrate) up to 9999 images. The summation mode not only reduces random noise in the final image by the square root of n but also acts as a gain

factor up to 9999 for weak signals. The summed image can also be divided by powers of two as high as 256, giving a variable-gain factor. The summing mode is therefore a combination of gain and averaging and is particularly useful for very low signal levels. Both the averaging and summing modes operate at 13.4 frames per second. The new image creation rate is the number of frames processed divided by the system frame rate. For an averaging level of 2, therefore, new images appear at 6.7 Hz, and summing eight frames takes 0.60 s.

COMMUNICATIONS

The light scanning microscope has a SCSI port to save images to any kind of hard drive or removable drive. The file format is TIFF, a universal standard compatible with most image processing programs. An ethernet port is included for connection to networks for storage, archiving, and communication across the Internet or other telecommunication systems.

Files that have been created or modified in TIFF format by image processing or desktop publishing programs can be read back into COSMIC and shown on the display for teaching or conferences. A "slide show" mode can sequentially display all the images on the disk, so this means that text slides or photographs can be displayed along with the microscope images.

OVERSAMPLING AND SUPERRESOLUTION

As mentioned earlier, oversampling is a process of collecting more electronic data points than actually exist in the optical image. If there are 1000 optical resolution elements in a scan line across the image, and that line is digitized into 3000 data points, then that is oversampling. The number of optical resolution elements in a field of view can easily be calculated by using the classical Airy disk formula

$$R_o = \frac{0.61 \lambda}{NA}$$

where λ is the wavelength, typically chosen as 550 nm, and NA is the numerical aperture of the objective. According to the Rayleigh resolution criterion, R_o is the minimum distance resolvable by a lens with a given numerical aperture. It is also the radius of the Airy disk, which is the smallest dot that can be imaged by an objective, and will be referred to as the Airy disk size or the optical resolution element size.

Figure 13.7 shows a 40x objective with a numerical aperture of 0.75 which has an Airy disk size of 0.447 μm that defines its resolution. In a normal field of view of diameter 500 μm , there are approximately 1100 optical resolution elements across the center of the field. The COSMIC CRT raster scans a 3:4 aspect ratio rectangle whose diagonal is 70% of the 500- μm circle diameter, resulting in a scanned area that is approximately 280 \times 210 μm , or about 626 \times 470 Airy disks.




Fig. 13.7. A 40x objective with a NA of 0.75 has an Airy disk size of $0.447 \mu\text{m}$ that defines its resolution. In a normal field of view of diameter $500 \mu\text{m}$, there are approximately 1100 optical resolution elements across the center of the field. The CRT raster scans a 3 : 4 aspect ratio rectangle whose diagonal is 70% of the 500- μm circle diameter, resulting in a scanned area that is approximately $280 \times 210 \mu\text{m}$, or about 626×470 Airy disks.

In order to determine the resolving power of the imaging system, it is necessary to take into account the sampling theory of Nyquist, who derived the number of samples needed to define a single point of a data curve. In terms of optics, to represent the information in a single optical resolution element (Airy disk), one must take two samples across the element. This is not really *oversampling* but *adequate sampling*, since 2 : 1 is the ratio needed for accurate sampling. Oversampling would therefore be anything over the 2 : 1 ratio. Examining Figure 13.7 again, it can be seen that for the 40x objective (with no zoom) 1280 points are digitized across the 626 Airy disks. The sampling ratio is right at 2:1, and it is reasonable to say that the data from the 40x objective are “just resolved.”

When the 40x field of view is zoomed by 3x, the number of optical resolution elements in the scanned area is reduced to 218×156 , but the data are still digitized into 1280×1024 data points. This is almost a 6 : 1 oversampling (Plasek et al., 1998). In reality, oversampling does give an improved image over “adequate sampling”; that is, the modulation transfer function improves up to a certain point, after which there is no improvement. The image detail increases with magnification up to a point, after which the image gets larger but no more detail is seen (“empty magnification”).

Each objective power has a different number of optical resolution elements within the field of view. Interestingly, lower power objectives have more elements across the field of view, posing a greater imaging challenge than high-power objectives. For example, a 2.5x objective has a scanned area of $4.300 \times 3.225 \text{ mm}$ that contains




Fig. 13.8. Superresolution effect is demonstrated with an image of the diatom *Pleurasigma angulatum*. The holes in the diatoms are approximately $0.25\ \mu\text{m}$ in diameter and hexagonal in shape. The full-color images were taken in white light. Obtaining such high-resolution images is very difficult with most other systems. Figure 13.8 (A) was taken with a 40x (0.75NA) objective. Figure 13.8 (B) was taken with a 100x (1.3NA oil) objective. However, the sample is in air under the cover slip, so the maximum usable NA for the 100x objective is less than 1.0. The condenser had an NA of 1.4 (oil). See color plates.

961 \times 721 Airy disks. The 1280 \times 1024 electronic data points are less than the 2 : 1 ratio needed to fully resolve the smallest element in the image. However, if the image is zoomed by 1.6, the scanned area covers 600 \times 400 Airy disks and all points in the field are fully resolved. The zoom control provides a convenient mechanism to trade field of view for full resolution whenever desired. Because the zoom factor will carry any objective into true empty magnification, the system resolution is limited by the objective's ability to create an Airy disk, not the ability to take the data from it. Therefore, it is reasonable to say that the limit of resolution of this digital microscope is the optics, not the electronics.

This color scanning microscope can also exhibit superresolution, where image detail far smaller than the calculated limit has been observed. The technique uses the oversampling feature coupled with COSMIC's summing, averaging, and background correction functions. Diatom images taken with a 40x objective (0.75 NA) have shown features as small as $0.05\ \mu\text{m}$, well below the Airy disk size of $0.447\ \mu\text{m}$. This superresolution effect is demonstrated in Figure 13.8, where an image is shown of *Pleurasigma angulatum*. The holes in these diatoms are approximately $0.25\ \mu\text{m}$ in diameter and are hexagonal in shape.

BRIGHTNESS AND COLOR CALIBRATION

The light transmitted through different objectives varies with the power of the objective. Significantly more light is transmitted through the 2.5x objective than through the 100x objective. In addition, the color balance of each objective is different, so the variation in red/green/blue balance must also be corrected. Calibration of brightness and

color balance to standard values is done by adjusting the gain setting of each PMT detector for each objective. A sensor embedded in the objective nose piece automatically selects the proper setup conditions when an objective is selected. These presets can be changed by the user if desired for special effects or if a new objective is placed on the instrument. The calibration procedure establishes a standard color balance and brightness for white, ensuring accurate color reproduction for stained samples.

COMPARISON WITH CAMERA SYSTEMS

Matching a CCD camera detector to a microscope image to ensure maximum resolution requires a different design procedure from COSMIC. CCD camera detectors consist of an x - y array of tiny photodiodes. Each photodiode represents one element of the picture and is called a pixel. Light striking the photodiode surface creates electrons that are read out as the signal. Since the pixels on a CCD detector are fixed, the size and spacing of the pixels must be compared with the Airy disk size to determine performance. For example, to evaluate a camera with a 12- μm square pixel, a simple calculation of the size of the Airy disk image on the detector will indicate resolution. A 100x objective (1.3 NA) has an Airy disk size of 0.256 μm in the sample plane. At the detector, the image of the Airy disk is $100 \times 0.256 \mu\text{m}$, or 25.6 μm . Since the Airy disk covers two 12- μm electronic pixels, the image will be fully resolved. If the objective is a 40x (0.85 NA), with an Airy disk of 0.394 μm , the image of the Airy disk is $40 \times 0.394 \mu\text{m}$, or 15.8 μm , which is not fully resolved by the detector, since the sampling ratio is only 1.31 : 1. On a CCD detector, there is no zoom feature available to trade field of view for resolution or to create oversampling. The situation is even worse for a 2.5x objective (0.75 NA), where the Airy disk size is 4.47 μm at the sample plane and 11.18 μm at the detector plane. The sampling ratio is only 0.93 : 1 for this example. To achieve good resolution, CCD camera systems must carefully take into account the pixel size on the CCD detector as well as the Airy disk size of each objective. The number of pixels does not determine resolution: it defines the field size. It is the size of the pixel compared to the size of the Airy disk on the image plane that determines resolution. Thus a camera with 3500×2700 pixels cannot claim a higher resolution than a camera with 1500×1200 unless each individual pixel is smaller.

The issues of field size and pixel size are very important in understanding the differences between the point scanning technique of COSMIC and CCD detectors. As with all semiconductor technologies, the size of CCD arrays has become smaller with improved manufacturing techniques. Sizes as commonly designated include 1 in. (12.8×9.6 mm), $\frac{2}{3}$ in. (8.8×6.6 mm), $\frac{1}{2}$ in. (6.4×4.8 mm), $\frac{1}{3}$ in. (4.8×3.6 mm), and $\frac{1}{4}$ in. (3.2×2.4 mm). The signal from a single photodiode pixel is directly proportional to the amount of light striking its surface. If the photodiode pixel is made smaller to get higher resolution (more pixels in the same area), then the signal decreases. A longer time is required to accumulate sufficient signal electrons, and the frame rate goes down accordingly. When light is divided between three CCD detectors to create color images, the loss of sensitivity is even greater. Frame rates for high-

resolution color CCD cameras are typically below 2 Hz, which is too slow for scanning around on a sample and makes focusing difficult. With these cameras it is necessary to focus visually through the microscope, then switch the light path to the camera and refocus and readjust the illumination. The process requires operating two instruments to get a digital image. The large image sizes of some cameras (viz. 3200×2400 pixels) creates an additional problem. The amount of data in one image can be so large that the time to transfer those data from the camera to a computer for display is significant, taking up to 1 min in some cases. Strategies to increase the frame rate for scanning and focusing include transmitting only a small fraction of the total image. The resulting tiny, low-resolution black-and-white images are not ideal for scanning or focusing. Therefore, the time required to establish the best focus, color levels, and intensity levels and store the image can be very long indeed.

In this scanning microscope the image has been standardized at 1280×1024 pixels, which means that the resolution is limited by the Airy disk size of the objective and not by the electronics or pixel size of a detector. With the scan rate of 13.4 Hz (frames per second), it is easy to focus and move the image. Indeed, one can obtain a large number of good images in a short time using this system—a major advantage from a user perspective. There is no reason why a larger image size could not be used, with the caveat that higher file size and slower image frequency would result. Of course, these factors are frequently based upon current computational capabilities, a factor that changes constantly.

There are several other significant differences between COSMIC and traditional camera-based systems. Most, if not all, of the following issues are not applicable to the white-light scanning microscope but may be significant problems for camera systems under certain operating conditions.

First, the regular array pattern of pixels in a CCD camera can cause the generation of false data. When the sample has a periodic structure whose spacing is about the same as that of the detector pixel spacing, a beat frequency can be generated from slight misalignments between the image and the detectors. The result is wavelike patterns appearing in the picture. These patterns, called aliases, or Moiré patterns, are not real in the image but are caused by the interaction of the image with the detector. These aliasing patterns can sometimes be eliminated by rotating the image against the detector. Aliasing can be a significant problem for certain kinds of images, particularly for automated image analysis. Aliasing problems can be greater for color cameras that have three CCD detectors, one each for the red, green, and blue portions of the light spectrum. The three detectors must be carefully aligned to give the same size image for each color and to prevent aliasing between the detectors.

Second, if a portion of a CCD detector is saturated with light, the excess electrons or the excess light may cause a signal from adjacent pixels. The visual effect may look like a star pattern with perpendicular arms and is called blooming. Electronics to detect and minimize blooming are usually built into the camera; however, the result is that the operation of the camera is temporarily halted.

Third, in an ordinary microscope, all points in the image field are illuminated simultaneously, which is called full-field illumination. In this scanning microscope, only a single point is illuminated at a time since illumination is provided by a flying

spot. When all the pixels are simultaneously illuminated, scattered light can reduce contrast from pixel to pixel and contribute to blooming.

The fourth issue is interpolation. Medical systems require better resolution than can be obtained from inexpensive mass-produced television-format CCD cameras. Some camera manufacturers have chosen to create artificial resolution by calculation. Several versions have been implemented. One method is to increase resolution simply by interpolation between pixels. In another case, one of the three color detectors is displaced one-half pixel horizontally and one-half pixel vertically. The color data for that physical location are computed from the nearest neighbors. In this technique, there are no physical locations with true information in any part of the image since all data must be computed. It is not clear how interpolation impacts the use of these detector systems, but the issue will have to be addressed eventually. Finally, because the sensitivity of PMTs is at least 1000 times greater than typical CCD detectors, the illumination intensity of the white-light scanning microscope can be 1000 times less than traditional light sources, reducing bleaching or photodamage to the specimen, an issue frequently overlooked in biological imaging.

OPERATIONAL CONSIDERATIONS

The most significant problems to be dealt with in a CRT-based spot-scanning instrument are the cost of the instrument and the aging characteristics of the CRT, which will decrease in light output over time. CRT lifetimes in the 10,000-h range are common, and replacement is not difficult or expensive, particularly when compared with most lasers used in microscopy. However, extended use of the zoom feature can cause preferential aging of the center of the CRT screen. Therefore, zoom time is limited by a timer circuit that disables the zoom after a preset interval (2 to 10 min.), thereby preventing the zoom feature from aging the center of the screen. The brightness and color balance calibration procedure continually corrects for CRT aging until the overall light level is unacceptable and a replacement is needed. Changing or adding different objectives requires calibration for that objective, but once performed, the calibration becomes an instrument preset.

A troublesome issue for those who must deal with microscope and camera systems is the difficulty of obtaining rapid and reproducible images from camera/computer systems. Combining the camera, imaging boards, and computer software is, in itself, an interesting challenge. In this system there are fewer variables and therefore a more consistent and stable image platform.

The clearest departure from traditional microscopy in this scanning microscope is the lack of a viewing ocular. This is a radical difference for microscopists who may feel uncomfortable not being able to view the specimen, except electronically. However, the format is most likely consistent with the future of digital imaging applications. Clearly, the whole concept of telepathology or remote evaluation of images implies that at least one participant will not have access to a traditional viewing ocular. Thus, the complete removal of the ocular may actually be advantageous for the originator of the images to know exactly what the remote viewer is seeing.

The COSMIC instrument is clearly the first fully integrated digital microscope system to offer real-time, full-color, high-resolution digital images. The unique features of zoom magnification, contrast enhancement and suppression, averaging and summing, background correction, and superresolution capabilities are major advantages. The lack of moving parts and alignment protocols is a powerful incentive for use in service roles such as pathology, where a robust, easy-to-operate, high-resolution instrument is desired.

FUTURE APPLICATIONS

Because of the nature of the optical pathway, transmitted fluorescence detection is likely to be a viable option using this technology. Currently, no commercial microscopes use transmitted fluorescence; they are all epifluorescent systems. Because bright-field and fluorescence in a simple, easy-to-use system is a desirable combination, addition of fluorescence to this scanning microscope will significantly enhance its potential in pathology.

REFERENCES

- Bartels PH, Layton J, Shoemaker RL (1984): Digital microscopy. *Monogr Clin Cytol* 9:28–61.
- Bibbo M, Bartels PH, Dytch HE, Puls JH, Pishotta FT, Wied GL (1983): High-resolution color video cytophotometry. *Cell Biophys* 5:61–69.
- Ingram M, Preston K, Jr (1970): Automatic analysis of blood cells. *Sci Am* 223:72–82.
- Ingram M, Preston K, Jr (1964): Importance of automatic pattern recognition techniques in the early detection of altered blood cell production. *Ann NY Acad Sci* 113:1066–1066.
- Mellors RC, Keane JF Jr, Papanicolaou GN (1952): Nucleic acid content of the squamous cancer cell. *Science* 116:265–265.
- Nagata H, Mizushima H (1998): World wide microscope: new concept of internet telepathology microscope and implementation of the prototype. *Medinfo 9 Pt 1*:286–289.
- Plasek J, Reischig J (1998): Transmitted-light microscopy for biology: A physicist's point of view. *Proc RMS* 33:196–205.
- Ploem JS (1988): Modern image analysis methods in hematology. *Nouv Rev Fr Hematol* 30:45–49.
- Ploem JS, van Driel-Kulker AMJ, Goyarts-Veldstra L, Ploem-Zaaijer JJ, Verwoerd NP, van der Zwan M (1986): Image analysis combined with quantitative cytochemistry. Results and instrumental developments for cancer diagnosis. *Histochemistry* 84:549–555.
- Prewitt JM, Mendelsohn ML (1996): The analysis of cell images. *Ann NY Acad Sci* 128:1035–1053.
- Singson RP, Natarajan S, Greenson JK, Marchevsky AM (1999): Virtual microscopy and the Internet as telepathology consultation tools. A study of gastrointestinal biopsy specimens. *Am J Clin Pathol* 111:792–795.
- Szymas J, Wolf G (1998): Telepathology by the internet. *Adv Clin Path* 2:133–135.
- Vazir MH, Loane MA, Wootton R (1998): A pilot study of low-cost dynamic telepathology using the public telephone network. *J Telemed Telecare* 4:168–171.